Accepted Manuscript

Endoscopy International Open

A new through-the-scope clip with anchor prongs is safe and successful for a variety of endoscopic uses: a case series (with videos)

John J Guardiola, Douglas K. Rex, Christopher C Thompson, Jeffrey Mosko, Marvin Ryou, Joyce Peetermans, Matthew J Rousseau, Daniel von Renteln.

Affiliations below.

DOI: 10.1055/a-2330-9803

Please cite this article as: Guardiola J J, Rex D, Thompson C C et al. A new through-the-scope clip with anchor prongs is safe and successful for a variety of endoscopic uses: a case series (with videos). Endoscopy International Open 2024. doi: 10.1055/a-2330-9803

Conflict of Interest: John J. Guardiola: Education – Olympus Corporation, Boston Scientific

Douglas K. Rex: Consultant - Olympus Corporation, Boston Scientific, Braintree Laboratories, Norgine, Medtronic, Acacia Pharmaceuticals; Research Support - Olympus Corporation, Medivators, Erbe USA Inc, Braintree Laboratories; Shareholder - Satisfai Health

Christopher C. Thompson: Apollo Endosurgery – Consultant/Research Support (Consulting fees/Institutional Research Grants), Bariendo – Founder/Board Member/Ownership Interest, BlueFlame Healthcare Venture Fund – Founder/General Partner, Boston Scientific – Consultant (Consulting fees)/Research Support (Institutional Research Grant), Medtronic – Consultant (Consulting Fees), ELLES – Founder/Board Member/Ownership Interest, Endoquest Robotics – Consultant, Institutional Research Grant, Enterasense Ltd – Founder, Consultant, Board Member, Ownership Interest, EnVision Endoscopy – Founder, Board Member, Consultant, Ownership Interest, ERBE – Institutional Research Grant, Fractyl – Consultant/Advisory Board Member (Consulting Fees)/Research Support, FujiFilm – Consultant/ Institutional Research Grant, GI Dynamics – Consultant (Consulting Fees)/ Research Support (Institutional Research Grant), GI Windows – Founder, Board Member, Ownership interest, Lumendi – Consultant/Institutional Research Grant, Olympus/Spiration – Consultant (Consulting Fees)/Research Support (Equipment Loans), Society for Metabolic and Bariatric Endoscopy (SMBE), Inc – Founder/President/ Ownership Interest, Softac – Consultant/Ownership Interest, USGI Medical – Consultant (Consulting Fees)/Advisory Board Member (Consulting fees)/Research Support (Institutional Research Grant), Xenter – Consultant/SAB/Ownership Interest

Jeffrey D. Mosko: research funding from Pendopharm and the Canadian Association of Gastroenterology and has received consultant or speaker fees from Boston Scientific Inc., Pendopharm, Vantage Endoscopy, Pentax, Fujifilm and Medtronic.

Marvin Ryou: Consultant and study support from Boston Scientific, Cook Medical, Olympus. Consultant for Fuji. Equity and Consultant for EnteraSense and GI Windows Surgical.

Joyce A. Peetermans: full-time employee of Boston Scientific

Matthew J. Rousseau: full-time employee of Boston Scientific

Daniel von Renteln: research funding from Boston Scientific Inc., ERBE Elektromedizin GmbH, Ventage, Pendopharm, Fujifilm and Pentax, and has received consultant or speaker fees from Boston Scientific Inc., ERBE Elektromedizin GmbH, Fujifilm and Pendopharm.

This study was supported by Boston Scientific Corporation (http://dx.doi.org/10.13039/100008497), Fonds de Recherche du Québec Santé Career Award

This is a PDF file of an unedited manuscript that has been accepted for publication. As a service to our customers we are providing this early version of the manuscript. The manuscript will undergo copyediting, typesetting, and review of the resulting proof before it is published in its final form. Please note that during the production process errors may be discovered which could affect the content, and all legal disclaimers that apply to the journal pertain.



Trial registration: NCT05653843, ClinicalTrials.gov (http://www.clinicaltrials.gov/), Prospective, Multi-Center

Abstract:

Background & Aims: Endoscopic through-the-scope clips (TTSC) are used for hemostasis and closure. We documented the performance of a new TTSC with anchor prongs.

Methods: We conducted a prospective case series of the new TTSC in 50 patients with an indication for endoscopic clipping at 3 hospitals in the USA and Canada. Patients were followed for 30 days after index procedure. Outcomes included defect closure and rate of serious adverse events (SAEs) related to the device or procedure.

Results: Fifty patients had 56 clipping procedures. Thirty-four procedures were clipping after EMR in colon (33) or stomach (1), 16 after polypectomy, 2 for hemostasis of active bleeding, and 1 each for fistula closure, post-POEM mucosal closure, or anchoring a feeding tube. Complete defect closure was achieved in 32/33 colon EMR defects and 21/22 other defects. All clips were placed per labeled directions for use. In 41 (82.0%) patients, prophylaxis of delayed bleeding was reported as an indication for endoscopic clipping. There were 3 instances of delayed bleeding. There were no device-related SAEs. The only technical difficulty was one instance of premature clip deployment.

Conclusion: A novel TTSC with anchor prongs showed success in a range of defect closures, an acceptable safety profile, and low incidence of technical difficulties. (ClinicalTrials.gov number, NCT05653843)

Corresponding Author:

Dr. John J Guardiola, Indiana University School of Medicine, Divsion of Gastroenterology and Hepatology, Indianapolis, United States, joguar@iu.edu

Affiliations:

John J Guardiola, Indiana University School of Medicine, Divsion of Gastroenterology and Hepatology, Indianapolis, United States Douglas K. Rex, Indiana University School of Medicine, Division of Gastroenterology and Hepatology, Indianapolis, United States Christopher C Thompson, Brigham and Women's Hospital, Division of Gastroenterology, Hepatology and Endoscopy, Boston, United States

[...]

Daniel von Renteln, University of Montreal, Gastroenterology, Montreal, Canada

This is a PDF file of an unedited manuscript that has been accepted for publication. As a service to our customers we are providing this early version of the manuscript. The manuscript will undergo copyediting, typesetting, and review of the resulting proof before it is published in its final form. Please note that during the production process errors may be discovered which could affect the content, and all legal disclaimers that apply to the journal pertain.



Abbreviations and acronyms

AGA: American Gastroenterological Association
APC: Argon Plasma Coagulation
COPD: chronic obstructive pulmonary disease
DAT: dual action tissue
EMR: endoscopic mucosal resection
FDA: US Food and Drug Administration
G-LST: granular lateral spreading tumor
IQR: interquartile range
NSAID: non-steroidal anti-inflammatory drug
OTSC: over-the-scope clip
POEM: peroral endoscopic myotomy
SAE: serious adverse event
SAS: Statistical Analysis System
SD: standard deviation
TTSC: through-the-scope endoscopic clip

INTRODUCTION

Through-the-scope endoscopic clips (TTSC) are used for a broad range of indications, including closure of gastrointestinal defects, with decreasing closure efficacy for larger defects.[1] Many types of TTSC are available, and physicians often select clips based on physical characteristics, cost or functional characteristics.[2] Defect characteristics and TTSC capabilities should be considered when choosing the closure device and technique.[1] Defect size is an important determinant of clip choice, i.e. small perforations that occur during endoscopic polyp resection can successfully be closed using TTSCs with favorable outcomes, while larger perforations are more expediently and effectively closed with large-caliber over-the-scope clips (OTSC) mounted on the outside of a clear cap affixed to the tip of the endoscope.[3]

A 2019 bench study tested 5 models of marketed TTSC with respect to rotatability, overshoot/whip, open/close precision, and tensile/closure strength in 4 different endoscope configurations: (1) straight, (2) duodenal sweep, (3) full retroflexion, and (4) across the duodenoscope elevator.[2] The authors mentioned that clip use in confined spaces or with thin tissue (e.g., Zenker's diverticulotomy, endoscopic submucosal dissection-induced defects or bleeding) probably calls for a smaller clip with precise open/close, while therapies in full retroflexion or across the duodenoscope elevator probably require precise clip functioning in the most strained configuration.[2] Perforation closure may require clips with enhanced ability to laterally manipulate tissue to appose defect edges, and large ulcers or fibrous ulcers may require a clip with higher compression force.[2] Because increased gripping ability might expand potential indications, we are testing a TTSC with anchor prongs. Use of this clip has been reported for closure of complex polyp resection defects in 3 patients.[4] For preliminary evidence of the functional capacity of the new TTSC, we present an analysis of our first 50 cases.

METHODS

Study Design

This was a case series of 50 consecutive patients with indication for endoscopic clipping, who received treatment with MANTIS clip (Boston Scientific Corporation, Marlborough, Massachusetts, USA). This clip was cleared by the US Food and Drug Administration (FDA) in August 2022.[5,6] The study was conducted at two USA sites and one Canadian site.

Patient population

Patients of any age scheduled for an indicated endoscopic clipping per local standard of practice were eligible to enroll. Exclusion criteria were enrollment in another study that would directly interfere with the current study or investigator's assessment that the subject was at risk for study device or procedure related complications. All centers obtained approval from their respective local ethics committees. All patients provided signed informed consent before nonemergent procedures. In emergent cases (e.g. perforation or acute bleeding), consent was performed after the procedure, but before data collection because preprocedural consent was unfeasible. Periprocedural management of antiplatelet and anticoagulant agents was per endoscopist discretion.

Study Visits

Index procedure and postprocedural follow-up

Baseline screening visit included informed consent, eligibility assessment, age, gender and relevant medical history. Intervention was placement of new endoscopic clip in gastrointestinal tract. During index procedure, the type of clipping procedure, number of clips placed, and lesion size was recorded. After index procedure, participants continued to receive medical care per standard of practice. Adverse events and device events were followed for 30 days. The last study visit was telephone interview at 30-day (\pm 5 days). For participants who did not complete this interview, the reason was recorded.

Outcomes

We evaluated: 1) successful closure of defects defined as no submucosa visible and clips < 1 cm apart, 2) delayed bleeding rate, 3) reinterventions, 4) rate of serious adverse events (SAEs) related to MANTIS clip or endoscopic portion of procedure, and 5) technical difficulties.

Statistical Analysis

Baseline characteristics, medical history, outcome measures, and adverse events were summarized using mean, median, standard deviation (SD), range for continuous variables (e.g. age, procedure times), and proportions for categorical variables. All analyses performed in SAS version 9.4.

RESULTS

Patient and procedure characteristics

Among 63 patients screened for study eligibility, 10 were excluded for not having an indication for clip placement, and 3 were excluded because the investigator deemed them at risk for study device or procedure-related complications per the Instructions for Use. Among the latter three patients, one defect was deemed too fibrotic, there was concern placing both jaws in the center of a hot EMR defect in the rectum of second, and there was concern placing both jaws of the MANTIS in submucosa to manage a vessel with ongoing bleeding in EMR defect in third patient. Standard TTSC were used in these cases. Among enrolled participants, mean age was 62.8 (range 28.0–80.0) years, most participants were male (n = 32, 64.0%) and white (n = 44, 88.0%). Sixteen (32.0%) were taking nonsteroidal anti-inflammatory drugs including aspirin, 6 (12.0%) were taking anticoagulants and 2 (4.0%) were taking antiplatelets at baseline (Table 1).

Types of procedures performed

Post-EMR prophylactic hemostasis was reason for clipping in 34 defects. EMR was performed for 33 colonic lesions and 1 gastric lesion. Some patients had multiple lesions and thus underwent EMR of multiple lesions or a combination EMR/hot snare polypectomy. 16 colon polyps were removed from 12 patients via polypectomy without EMR – 15 using hot snare and 1 using cold snare. These resection defects were then clipped. Clipping was performed in two patients (4.0%) for hemostasis of active bleeding. Clipping was performed in one patient (2.0%) each for gastrocutaneous fistula closure[7], post peroral endoscopy myotomy (POEM) closure of the mucosal incision, or anchoring of a feeding tube (**Table 2**).

Study Outcomes Successful Defect Closure

Successful defect closure was completed in 53/55 (96.4%) lesions in 47 total patients in this study. The case in which endoscopic clipping was used for successful feeding tube anchoring is excluded here. Colonic post-EMR defects (33) were most common lesion. All clips were placed per labeled directions for use. A post-EMR defect measuring 30 mm in maximum diameter x 5 mm in minimum diameter in the cecum was unable to be completely closed due to difficult positioning. A rectal polypectomy defect measuring 20 x 20 mm was incompletely closed due to inability to appose tissue.

Among the 55 defects, median maximum defect size was 18.0 mm (range 1.2-60.0 mm), and median minimum defect size was 10.0 mm (range 1.0-40.0 mm). The mean number of study clips used was 2.6 ± 1.8 (range 1.0-9.0), and the mean number of other (non-study) clips used was 1.0 ± 1.9 (range 0-9.0) per procedure. The mean total number of clips used during a procedure was 3.6 ± 2.6 (range 1.0-11.0).

For colon post-EMR defects, 32/33 defects were completely closed. The median minimum defect diameter was 10 mm (range 5-40 mm) and the median maximum defect diameter was 25 mm (range 8-45 mm) for these defects. For post-EMR defects of colon, the mean number of clips placed was 4.4±2.6 (range 1 to 11, including an average of 3.0 study clips and 1.3 non-study clips). Due to some defects lying over folds and other factors, the defect size was smaller than the actual lesion removed in some cases.

For non-colon EMR defects (excluding using study clip to affix jejunal tube), 21/22 were successfully completely closed. The median minimum diameter was 8 mm (range 1-40 mm) and the median maximum diameter was 10 mm (range 1.2-60 mm). For all non-EMR of colon defects, the mean number of clips placed was 2.7±2.4 (range 1 to 11, including an average of 2.1 study clips and 0.6 non-study clips).

Delayed Bleeding Prophylaxis

In 41 (82.0%) patient cases, prophylaxis to reduce risk of delayed bleeding was reported as indication for endoscopic clipping. Delayed bleeding from clipped intervention sites occurred in 3 of 56 (5.4%) total interventions in 50 patients in the study. For post-EMR colon defects specifically, delayed bleeding occurred in 1/33 (3.0%) clipped defects. The one case of delayed bleeding was previously mentioned incompletely closed cecal defect. For all other defects, delayed bleeding occurred in 2/22 (9.0%) clipped defects. These cases were of the previously mentioned incompletely closed rectal polypectomy defect and the post-POEM defect.

Serious Adverse Events Related to the MANTIS Clip

Of the 50 enrolled patients, none of the patients had AEs or SAEs related to the MANTIS clip nor any other clip.

Study Clip Technical Success

There was one instance of the clip deploying prematurely. There were zero instances of the clip failing to deploy, failing to release from the catheter, being difficult to release from the catheter, failing to open or close, being unable to rotate, failing to advance from the sheath, being difficult to deploy, having the clip arms or anchor prongs bent, having poor bite of tissue, or any other events that could be seen as technical failure of the clip itself.

Serious Adverse Events Related to the Endoscopic Procedure and Reinterventions Five patients had a total of 4 SAEs and 1 nonserious AE related to the endoscopic portion of the procedure within a median of 12 (range 0 - 14) days after the index procedure (Table 3). The first patient had POEM procedure with a related SAE and 2 reinterventions. The patient was found to have esophageal leak on imaging on day 0 but site was contained so patient was treated with antibiotics and nothing by mouth. This patient had two repeat upper endoscopies. Patient had bleeding from esophageal mucosectomy site and an esophageal ulcer. The first repeat endoscopy performed for bleeding on day 13 after procedure showed a clot near the clipping site; the study clip was removed for better visualization and the site was treated with a hemostatic spray and placement of a non-study clip. Two days later, this patient had another procedure to treat an esophageal ulceration using hemostatic spray. No blood transfusion was required. The second patient had a related SAE and reintervention due to bleeding on day 12 after EMR of cecal lesion that was incompletely closed due to difficult access. This patient underwent colonoscopy. The endoscopist surveilled the affected site, but no treatment was needed. No blood transfusion was required. A third patient who had a rectal hot snare polypectomy site that was not completely closed also had a bleed on post procedure day 6 that resolved (nonserious AE with a reintervention). Subsequently, colonoscopy was performed and no intervention was needed. No blood transfusion was needed. Notably, a fourth patient had an unrelated SAE with a reintervention after melena occurred on day 1 after procedure. An upper endoscopy demonstrated bleeding not related to the previous endoscopic procedure nor the study clip. This bleeding was treated with a non-study clip. A fifth patient had a related SAE of postpolypectomy electrocoagulation syndrome that did not require reintervention.

Descriptions of how clip was used

The "grasp-and-drag" technique for a standard EMR defect

After hot snare piecemeal EMR of a lateral spreading lesion in the colon, several non-MANTIS clips had already been placed. Video 1 demonstrates the MANTIS clip was opened, and normal mucosa on one edge of the lesion was grasped. With the MANTIS clip closed, the colonoscope was used to lift the MANTIS clip with grasped tissue to the opposite side of the lesion. The MANTIS clip was then opened, with the anchor prongs maintaining the grasp of tissue. With the jaws open, the MANTIS clip was gently pushed to the mucosa to maximize use of its jaws, then closed approximating defect edges, then deployed. (Video – Part 1)

The "open-jaw" technique for a scarred defect

The lesion was a 50-mm granular lateral spreading tumor (G-LST) in the ascending colon that had undergone a previous attempt at resection before referral to our center. The lesion was assessed and resected with piecemeal hot EMR. Evidence of the previous resection attempt was seen at the haustral folds as scarred-down tissue extending across the defect. Hot avulsion was used to remove flat and fibrotic polyp tissue. Snare tip soft coagulation was performed to the edges of the defect.

The MANTIS clip's anchor prongs being used to drag tissue to approximate the defect without closing and reopening the jaws ("open jaw technique"), or slipping from the tissue as standard through the scope clips might when used to close a wide defect. This technique is useful when the defect is approached tangentially. Here, the MANTIS clip was rotated so that when opened, the orientation was vertical with the top jaw between 10 and 2 o'clock. The bottom anchor prongs were used to grasp normal tissue at the distal (anal side) edge of the lesion, and the colonoscope was advanced to lift and push the mucosa on the anal side of the defect toward the mucosal edge on the defect's cecal side. The top anchor prongs were then used to grasp normal tissue at the proximal edge of the lesion. The clip was pushed gently to "bury" the clip and use the full length of the jaws before deployment. This sequence was repeated. Defect was closed with total 6 MANTIS clips and 1 standard TTSC. The MANTIS clip was easily maneuvered for placement between previously deployed clips. Alternatively, fewer MANTIS clips can be used to bring the defect edges closer, followed by closure of the remaining defect using standard TTSC. (Video – Part 2)

The "open-jaw" technique for a broad-based defect

This shows a large submucosal defect after piecemeal hot EMR of a 45-mm G-LST in the cecum. The "open-jaw" technique was again utilized to grasp normal tissue on the anal side of the lesion and move it to the cecal edge of lesion. Three MANTIS clips were used followed by three standard TTSC to completely close this lesion. (Video – Part 3)

DISCUSSION

A new TTSC with anchor prongs had a high rate of success for defect closure, low rate of delayed bleeds, and high rate of technical success regarding clip deployment in this case series of 50 patients across 3 centers in the US and Canada in a variety of cases. There were zero adverse events related to the clip itself and it was technically successful in all but one instance because of premature deployment. The low rate of related SAEs was consistent with published safety data on endoscopic clips.[8] We demonstrate two techniques for successful use of this technique.

Due to the tissue apposition ability of the MANTIS clip, it may be able to reduce the number of standard TTSCs used during defect closure or allow for the closure of larger defects without OTSC or suturing devices. Recently, another novel TTSC has been introduced – the dual action tissue (DAT) clip (Micro-Tech Endoscopy, USA, Ann Arbor, Michigan, USA).[9] This clip features two independent arms and is efficacious in closing large resection defects. It cannot be rotated so requires manipulation of the scope dials for positioning and operating the independent jaws can be difficult with endoscope looping.[10] The MANTIS is operated similarly to a standard TTSC. It may have a lower learning curve for both endoscopist and technician than alternative devices. Comparative studies for endoscopic defect closure between these clips and other devices such as suturing are warranted.

Data sharing

The data, analytic methods, and study materials for this study may be made available to other

researchers in accordance with the Boston Scientific Data Sharing Policy

(http://www.bostonscientific.com/en-US/data-sharing-requests.html).

References

- 1. McCarty TR, Jirapinyo P. Endoscopic Closure: Tools and Techniques. Gastrointest Endosc Clin N Am 2023; 33: 169-182 DOI: 10.1016/j.giec.2022.08.003
- 2. Wang TJ, Aihara H, Thompson AC et al. Choosing the right through-the-scope clip: a rigorous comparison of rotatability, whip, open/close precision, and closure strength (with videos). Gastrointest Endosc 2019; 89: 77-86 e71 DOI: 10.1016/j.gie.2018.07.025
- Elmunzer BJ. Just Clip It: Endoscopic Clipping in the 21st Century. Am J Gastroenterol 2016; 111:
 6-8 DOI: 10.1038/ajg.2015.382
- 4. Wei MT, Friedland S. Use of anchor pronged clips to close complex polyp resection defects. VideoGIE 2023; 8: 245-246 DOI: 10.1016/j.vgie.2023.02.002
- 5. US Food and Drug Administration. 510(k) number K221124: Mantis Clip (August 4, 2022). https://www.accessdata.fda.gov/cdrh_docs/pdf22/K221124.pdf. Accessed August 15, 2023.:
- 6. Boston Scientific Corporation. Directions for use: MANTIS Clip. <u>https://www.bostonscientific.com/content/dam/bostonscientific/endo/portfolio-group/</u> <u>mantis/MANTIS Short DFU 16Sep2022.pdf</u>. Accessed on August 15, 2023.:
- Ego M, Abe S, Nonaka S et al. Endoscopic Closure Utilizing Endoloop and Endoclips After Gastric Endoscopic Submucosal Dissection for Patients on Antithrombotic Therapy. Dig Dis Sci 2021; 66: 2336-2344 DOI: 10.1007/s10620-020-06508-8
- 8. Ramai D, Deliwala SS, Mozell D et al. Adverse Events of Endoscopic Clip Placement: A MAUDE Database Analysis. J Clin Gastroenterol 2023: DOI: 10.1097/MCG.00000000001822
- 9. Yang D, Kadkhodayan K, Arain MA et al. Novel dual-action tissue through-the-scope clip for endoscopic closure. VideoGIE 2022; 7: 345-347 DOI: 10.1016/j.vgie.2022.07.004
- 10. Mohammed ASS, Robalino Gonzaga E, Hasan MK et al. Low Delayed Bleeding and High Complete Closure Rate of Mucosal Defects with the Novel Through-The-Scope Dual Action Tissue Clip After

Endoscopic Resection of Large Non-Pedunculated Colorectal Lesions. Gastrointest Endosc 2023: DOI: 10.1016/j.gie.2023.07.025

- 11. Motchum L, Levenick JM, Djinbachian R et al. EMR combined with hybrid argon plasma coagulation to prevent recurrence of large nonpedunculated colorectal polyps (with videos). Gastrointest Endosc 2022; 96: 840-848.e842 DOI: 10.1016/j.gie.2022.06.018
- 12. Djinbachian R, Pohl H, Rex DK et al. Thermal ablation after endoscopic mucosal resection of large colorectal polyps: not only the margins, but also the base? Gut 2023; 73: 12-15 DOI: 10.1136/gutjnl-2023-331057

reserved		
All rights		
yright. /		
copyri		
by		
protected		
cle is		
arti		
his		

Table 1. Baseline patient characteristics (n=50)

Characteristic	Mean ± SD (range) or n/N (%)	
Mean age, yrs	62.8±11.9 (28.0, 80.0)	
Gender: Male	32/50 (64.0%)	
Ethnicity		
Not Hispanic or Latino	45/50 (90.0%)	
Hispanic or Latino	1/50 (2.0%)	
Not disclosed	4/50 (8.0%)	
Race		
White (Caucasian)	44/50 (88.0%)	
Black or African American	3/50 (6.0%)	
American Indian or Alaska Native	1/50 (2.0%)	
Hispanic or Latino	1/50 (2.0%)	
Not disclosed	1/50 (2.0%)	
Medical History*		
No medical history conditions	36/50 (72.0%)	
Coronary artery disease	9/50 (18.0%)	
Bleeding risk	5/50 (10.0%)	
Chronic kidney disease	4/50 (8.0%)	
Congestive heart failure	4/50 (8.0%)	
Chronic obstructive pulmonary disease	1/50 (2.0%)	
Esophageal dysphagia	1/50 (2.0%)	
Known elevated bilirubin	1/50 (2.0%)	
Liver disease/end stage liver disease	1/50 (2.0%)	
Lung cancer	1/50 (2.0%)	
Recurrent pneumonia	1/50 (2.0%)	
Actively taking NSAIDs	16/50 (32.0%)	
Actively taking anticoagulants	6/50 (12.0%)	
Apixaban	3/50 (6.0%)	
Lovenox	2/50 (4.0%)	
Warfarin	1/50 (2.0%)	
Actively taking antiplatelets	2/50 (4.0%)	
(clopidogrel)	2/30 (4.070)	

POEM

peroral endoscopy myotomy

*Each patient had 1 or more of the listed medical conditions; rows are not mutually exclusive.

Table 2. Procedural details (N=50)

Characteristic	Mean ± SD or median (range) or % (n/N patients)
Reason for clipping procedure	
Post-endoscopic mucosal resection (EMR)*	66.0% (33/50)
Polypectomy	24.0% (12/50)
Hemostasis of active bleeding	4.0% (2/50)
Fistula closure	2.0% (1/50)
Post-POEM mucosal incision closure	2.0% (1/50)
Anchoring/affixing of jejunal feeding tube	2.0% (1/50)
Additional modalities used during procedure	
None	86.0% (43/50)
Argon Plasma Coagulation with or without ligation^	4.0% (2/50)
Endoscopic Suturing - X-tack	4.0% (2/50)
Snare tip soft coagulation for active bleeding	4.0% (2/50)
Hemostatic grasper	2.0% (1/50)
Median maximum lesion diameter (range), mm	18.0 (1.2, 60.0)
Median minimum lesion diameter (range), mm	10.0 (1.0,40.0)
Mean number of study clips used per procedure (range)	2.6±1.8 (1.0,9.0)
Mean number of non-study clips used per procedure (range)	1.0±1.9 (0.0,9.0)
Mean total number of clips used per procedure (range)	3.6±2.6 (1.0,11.0)

*In the post-EMR category, 1 patient had 3 lesions, 4 patients had 2 lesions, and 28 patients had 1 lesion (total 39 lesions). Patients in all other categories had 1 lesion/clipping procedure each.

^APC was used for margin and base ablation of a colonic EMR defect[11,12] and combined with ligature device for gastrocutaneous fistula closure[7]

Table 3. Serious adverse events related to the endoscopic portion of the procedure

	Number of SAEs	Number of patients (n/N) (%)
Any serious adverse event	5	3/50 (6.0%)
Bleeding	2	2/50 (4.0%)
Esophageal perforation	1	1/50 (2.0%)
Post-polypectomy electrocoagulation syndrome	1	1/50 (2.0%)
Ulceration	1	1/50 (2.0%)

VIDEO LEGEND

Video. Part 1: "Grasp and drag" technique to close a defect after hot snare piecemeal endoscopic mucosal resection of a lateral spreading tumor in the colon. Part 2: "Open jaw" technique to prophylactically close a scarred defect of a hot endoscopic mucosal resection of a 50-mm scarred granular lateral spreading tumor in the ascending colon. Part 3: Repeated "open jaw" technique to prophylactically close a broad-based defect after piecemeal endoscopic mucosal resection of a 45-mm granular lateral spreading tumor in the cecum