This document was downloaded for personal use only. Unauthorized distribution is strictly prohibited.

Cut Guide Wire During Endoscopic Sphincterotomy

E. Parlak, S. O. Kuran, E. Altiparmak, Department of Gastroenterology, Türkiye Yüksek İhtisas Hospital, Ankara, Turkey

Corresponding Author

E. Parlak, M.D. Dept. of Gastroenterology Türkiye Yüksek İhtisas Hospital 06100 Sihhiye Ankara Turkey

Fax: +90-312-312-4122 E-mail: Eparlak@ato.org.tr



Figure 1 A 67-year-old man with cancer of the pancreatic head was admitted to our endoscopic retrograde cholangiopancreatography (ERCP) unit for biliary decompression. The biliary system was canalized using a guide wire (0.021-inch Black Wire, Bard Angiomed, Germany), and endoscopic sphincterotomy was then performed over the same guide wire. At the end of the sphincterotomy, it was noted that the guide wire had been cut in the common bile duct



Figure 2 The tip of guide wire, which was in the common bile duct, was brought out through the papilla with the help of a stone extraction balloon catheter



Figure 3 The tip of the guide wire in the duodenum was grasped with a snare loop and removed



Figure 4 The tip of the guide wire in the duodenum was grasped with a snare loop and removed