SPOTLIGHT 1789

SYNLETT Spotlight 165

Dinitrogen Tetroxide: N₂O₄

Compiled by Morteza Shiri

This feature focuses on a reagent chosen by a postgraduate, highlighting the uses and preparation of the reagent in current research

Morteza Shiri was born in 1978 in Hamadan, Iran. He obtained his B.Sc. in Applied Chemistry from Bu-Ali Sina University, Hamadan and his M.Sc. in Organic Chemistry from Shiraz University under the guidance of Professor Habib Firouzabadi. Presently he is a PhD student at the Bu-Ali Sina University, Hamadan under the supervision of Professor Mohammad Ali Zolfigol. His research field is application of new solid-supported reagents in organic reactions and synthesis of new indole derivatives.

College of Chemistry, Bu-Ali Sina University, 6517838683 Hamadan, Iran E-mail: mshiri@basu.ac.ir



Introduction

For over a century, dinitrogen tetroxide has found wide application in organic synthesis, such as nitration of aromatic compounds, ¹ nitrosation of amines, ² preparation of thionitrite³ and sulfinyl nitrites, ⁴ oxidation of olefins ⁵ and dethioacetalization reactions. ⁶ This compound is commercially supplied at low price in a cylinder, and can be used directly or as liquid (bp 21 °C), collected by transfer distillation into an ice-cooled vessel. Apart from difficul-

ties in handling the poisonous and corrosive dinitrogen tetroxide, the biggest disadvantage of utilizing gaseous N_2O_4 in organic reactions is its high reactivity which usually causes undesired side reactions. In order to overcome the above mentioned limitations some reports are published on the use of N_2O_4 complexes of organic, polymeric and inorganic compounds as useful reagents in organic reactions.

Abstract

 N_2O_4 was supported on the cross-linked polyvinylpyrrolidone (PVP) by Iranpoor et al., which afforded a stable polymeric reagent. Thiols were converted to S-nitrosothiols (thionitrites) using this new nitrosating agent in n-hexane or CHCl $_3$ at 10 °C. With this reagent, thiols were also converted into their corresponding disulfides. Selective oxidation of sulfides to sulfoxides and disulfides to thiosulfonates can also be achieved by this reagent at room temperature. By using an excess of the reagent, the selective one-pot synthesis of thiosulfonates from thiols at room temperature was also performed. 8a

RSH
$$\xrightarrow{\text{FVF-N}_2\text{O}_4}$$
 RSNO $\xrightarrow{\text{r.t.}}$ RSSR

RSH $\xrightarrow{\text{PVP-N}_2\text{O}_4}$ excess, r.t. RSO₂SR

RSSR $\xrightarrow{\text{PVP-N}_2\text{O}_4}$ r.t. RSOR

RSOR RSOR

Zolfigol and co-workers reported a practical and efficient method for N-nitrosation of secondary amines using $[NO^+\cdot Crown \cdot H(NO_3)_2^-]$, an ionic complex of 18-crown-6 with N_2O_4 gas . This complex, which is easily prepared and handled, is relatively stable and an efficient source for the delivery of nitrosonium ion (NO^+) under mild and homogeneous conditions. 18-Crown-6 can be recycled and reused.

$$\begin{array}{c} \text{RR'NH} & \xrightarrow{ \left[\text{NO}^+ \cdot \text{crown} \cdot \text{H}(\text{NO}_3)_2^- \right] } \\ \hline & \text{CH}_2\text{Cl}_2, \text{ r.t.} \end{array}$$

Firouzabadi et al. provided a simple method for the direct, regioselective iodination and bromination of benzene, naphthalene and other activated aromatic compounds using molecular iodine and bromine and their sodium salts in the presence of the stable $\rm Fe(NO_3)_3 \cdot 1.5N_2O_4/charcoal~system.^{9a}$

$$R \stackrel{R'}{=} = \frac{Fe(NO_3)_{3}\cdot 1.5 \ N_2O_4/charcoal, \ X_2 \text{ or NaX}}{anhyd \ CH_2Cl_2, \ r.t.} \qquad R \stackrel{R'}{=} X$$

$$X = I, Br$$

$$R, R' = H, Ar, Alk$$

SYNLETT 2006, No. 11, pp 1789–1790 Advanced online publication: 04.07.2006 DOI: 10.1055/s-2006-944218; Art ID: V17005ST © Georg Thieme Verlag Stuttgart ⋅ New York 1790 SPOTLIGHT

The C-2 lithiation of *N*-Boc- and *N*-(phenylsulfonyl)indoles, followed by reaction with dinitrogen tetroxide at low temperature affords the corresponding 2-nitroindoles in 63–78% yields. Deprotection of the *N*-Boc-2-nitroindoles with trifluoroacetic acid gives 2-nitroindole and 3-methyl-2-nitroindole in essentially quantitative yields.¹⁰

R = H, Me

1)
$$f$$
-BuLi, -78 °C
2) N_2O_4 , -120 °C to r.t.

N N NO₂

Various aziridines were reacted with dinitrogen tetroxide in the presence of Et₃N in anhyd THF to give the corresponding ethylenes in good to excellent yields.¹¹

[60]Fullerene was selectively and efficiently nitrated with dinitrogen tetroxide to yield monodisperse hexanitro[60]fullerene. The use of this compound as reactive precursor in the synthesis of organo-amino derivatives of C_{60} was also demonstrated. The allylic *tert*-nitro moieties in hexanitro[60]fullerenes were found to be excellent leaving groups in nucleophilic substitutions by amino nucleophiles, such as anilines, leading to the formation of hexanilino[60]fullerenes.¹²

$$C_{60} \xrightarrow{N_2O_4} \left(NO_2\right)_4 + \left(NO_2\right)_6$$

$$\left(NO_2\right)_6$$

Reaction of N_2O_4 with 1,2-dimethylcyclohexene led to the corresponding dinitro compound with high *trans* stereoselectivity (>30:1). Catalytic hydrogenation of 1,2-dimethyl-1,2-dinitrocyclohexene with $Pd(OH)_2$ on carbon afforded a quantitative yield of (d,l)-1,2-diamino-1,2-dimethylcyclohexane, which was resolved to pure enantiomers with mandelic acid. This approach constitutes a simple and potentially general route to interesting new chiral auxiliaries.¹³

Carboxamides and sulfonamides were reacted with dinitrogen tetroxide to give the corresponding acids, while such as *N*-bromoamides and hydroxamic acids were also converted to the corresponding acids together with their anhydrides.¹⁴

$$\begin{array}{c} O \\ \parallel \\ R\text{-C-NHX} \\ \text{ or } \\ R\text{-SO}_2\text{NH}_2 \\ \end{array} \begin{array}{c} O \\ \parallel \\ \parallel \\ R\text{-CO}_2\text{H} + R\text{-C-O-C-R} \\ \text{ or } \\ R\text{-SO}_3\text{H} \\ \end{array}$$

A novel class of stable, mild and size-shape-selective nitrosating agents for secondary amides was introduced by Zyryanov et al. These are based on reversible entrapment and release of reactive nitrosonium species by calix[4]arenes. The NO⁺ encapsulation controls the reaction selectivity.¹⁵

References and Notes

- (1) Squadrito, G. L.; Fronczek, F. R.; Watkins, S. F.; Church, D. F.; Pryor, W. A. *J. Org. Chem.* **1990**, *55*, 4322.
- (2) Boyer, J. H.; Pillai, T. P. J. Chem. Soc., Perkin Trans. I 1985, 1661.
- (3) Oae, S.; Shinhama, K.; Kim, Y. H. J. Chem. Soc., Chem. Commun. 1977, 407.
- (4) Oae, S.; Shinhama, K.; Kim, Y. H. Tetrahedron Lett. 1979, 35, 3307.
- (5) Giamalva, D. H.; Kenion, G. B.; Church, D. F.; Pryor, W. A. J. Am. Chem. Soc. 1987, 109, 7059.
- (6) Mehta, G.; Uma, R. Tetrahedron Lett. 1996, 37, 1897.
- (7) Zolfigol, M. A.; Zebarjadian, M. H.; Chehardoli, G.; Keypour, H.; Salehzadeh, S.; Shamsipur, M. J. Org. Chem. 2001, 66, 3619.
- (8) (a) Iranpoor, N.; Firouzabadi, H.; Pourali, A. Tetrahedron 2002, 58, 5179. (b) Iranpoor, N.; Firouzabadi, H.; Pourali, A. Synthesis 2003, 1591.
- (9) (a) Firouzabadi, H.; Iranpoor, N.; Shiri, M. *Tetrahedron Lett.* **2003**, *44*, 8781. (b) Iranpoor, N.; Firouzabadi, H.; Pourali,

- A. Synlett 2004, 347. (c) Iranpoor, N.; Firouzabadi, H.; Pourali, A. Phosphorus, Sulfur Silicon Relat. Elem. 2006, 181, 473. (d) Iranpoor, N.; Firouzabadi, H.; Heydari, R.; Shiri, M. Synth. Commun. 2005, 35, 263. (e) Iranpoor, N.; Firouzabadi, H.; Pourali, A. Synth. Commun. 2005, 35, 1517. (f) Iranpoor, N.; Firouzabadi, H.; Pourali, A. Synth. Commun. 2005, 35, 1527. (g) Iranpoor, N.; Firouzabadi, H.; Heydari, R. Synth. Commun. 2003, 33, 703.
- (10) Jiang, J.; Gribble, G. W. Tetrahedron Lett. 2002, 43, 4115.
- (11) Lee, K.; Kim, Y. H. Synth. Commun. 1999, 29, 1241.
- (12) Anantharaj, V.; Bhonsle, J.; Canteenwala, T.; Chiang, L. Y. J. Chem. Soc., Perkin Trans. 1 1999, 31.
- (13) Zhang, W.; Jacobsen, E. N. Tetrahedron Lett. 1991, 32, 1711.
- (14) Kim, Y. H.; Kim, K.; Park, Y. J. Tetrahedron Lett. 1990, 31, 3893.
- (15) (a) Zyryanov, G. V.; Rudkevich, D. M. Org. Lett. 2003, 5, 1253. (b) Zyryanov, G. V.; Rudkevich, D. M. J. Am. Chem. Soc. 2004, 126, 4264.